



# Antepartum Fetal Monitor BT-300



# Antepartum Fetal Monitor BT-300

- Compact and light
- Easy-to-use
- High sensitivity, probe-frequency 1 MHz
- Simultaneous external twin FHR with dual pulsed Doppler
- External uterine-contraction
- Fetal Movement
- Ergonomic designed event marker



## Specifications

### FHR

Dual Pulsed Doppler  
 Ultrasound Frequency: 1 MHz  
 Intensity:  $\leq 10 \text{ mW/cm}^2$   
 FHR Range: 50 - 240 bpm  
 FHR Accuracy:  $\pm 2\%$  of Range

### FM

Fetal Movement

### UC

External Type  
 Reference (zero) Control  
 Measurement Range: 0-99 units

### Printer

Thermal Array Type  
 Print Speed: 1, 2, 3 cm/min  
 Paper Feeding Function

### Display

7-Segment LED  
 3 Channels (FHR I, FHR II, UC)

### Sound

Doppler Sound  
 Alarm Sound  
 Information Sound

### Set-up

Alarm Upper / Lower Limit Value  
 Print Speed, Print Contrast

### Function

Mark Function  
 FHR II Offset Function  
 Auto Print Function

### Power

Input: AC 100-240 V (50/60 Hz)  
 Output: DC 16 V, 2.8 A  
 Consumption: 80 VA

### Dimension

191 (H) x 187 (L) x 201 (D) mm

### Standard Accessory

Ultrasound Doppler Probes (2)  
 UC Probe  
 Event Marker Jack  
 Print Paper  
 Power Adaptor & Cord  
 Ultrasound Gel  
 Probe Belt  
 Operation Manual

**Warranty** 2 Years